

GMP Flavour Standards

Product Information Sheet

Product Name: Geraniol Flavour Standard

Product Number: 190-0162-5-GMP

Document Code: PIS Geraniol 0003

Revision Number: 0003

Revision Date: 15th February 2016

Product Description

Geraniol GMP Flavour Standards are used in the sensory training and validation of sensory panels in the beer and beverage industry.

When placed into beer or another liquid this flavour standard imparts a rose water or rose hips like flavour. Sensory panellists should learn to taste and smell this flavour as part of a regular sensory quality control training program.

Geraniol GMP Flavour Standards are provided in powder form in single dose gelatine capsules within a pharmaceutical grade blister pack for improved product stability and ease of use.



Product Quality - GMP Flavour Standards

All FlavorActiV GMP Flavour Standards are produced to pharmaceutical grade quality control systems and standards within our dedicated production facility.

Our flavour centre uses Good Manufacturing Practices (GMP) and this quality system has been audited and certified by the Medicines and Healthcare Regulatory Authority (MHRA).

Our dedicated flavour centre is built to pharmaceutical grade specifications around the concept of handling flavours in liquid and powder forms in a controlled environment. This eliminates any flavour cross contamination and any extraneous taints and chemicals from entering our GMP Flavour Standards that may detract from their auglity.

This ensures that each FlavorActiV GMP Flavour Standard is produced to highest quality possible.

Product Use

- 1. To use a flavour standard push the gelatine capsule through the protective silver foil.
- 2. Prepare the volume of beverage you would like to impart the flavour to recommended 1 litre. Allow any foam to reduce.
- 3. Gently tap the gelatine capsule with your finger and remove the top half.
- 4. Pour slowly and evenly into the beverage.
- 5. Swirl the beverage for a few moments to distribute the flavour.
- 6. Pour samples and taste the beverage with your sensory panel





GMP Flavour Standards

Product Information Sheet

Product Name: Geraniol Flavour Standard

Product Number: 190-0162-5-GMP

Document Code: PIS Geraniol 0003

Revision Number: 0003

Revision Date: 15th February 2016

Technical Specifications:	
Product form:	Flavour standards in pharmaceutical grade blister packaging.
Flavour Imparted:	Geraniol
FEMA flavour number:	2507
EC chemical number:	203-377-1
CAS chemical number:	106-24-1
Amount of flavour per flavour standard:	1.2 mg (three times sensory flavour threshold)
Excipents per capsule:	Sucrose, B-cyclodextrin
Average volumetric contents per flavour standard:	500 mg
Solubility:	Soluble in all liquid beverages and foodstuffs
Quality control system:	Pharmaceutical GMP quality systems
Country of manufacture:	United Kingdom

Product Health Comments:	
Acceptable daily intake – beverage (ADI) -	FEMA Generally Regarded As Safe (GRAS) flavouring.
Microbiological contaminants - yeast and moulds	<100 cfu / capsule (<100 cfu per 100 ml in 1 litre)
Microbiological contaminants - bacteria	<100 cfu / capsule (<100 cfu per 100 ml in 1 litre)
Additional Comments:	N/A
Allergen Advice:	This GMP flavour standard is free from all allergens listed under EU FIC requirements and is produced in a pharmaceutical allergen free environment under strict controls with no risk of cross contamination.
	Those with a perfume allergy should not consume beverages prepared with this GMP Flavour Standard.

Other Comments:	
Disclaimer:	The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. FlavorActiV TM shall not be held liable for any damage resulting from handling or from contact with the above product.
Manufacturer Notice:	FlavorActiV™ flavour standards are manufactured by Quay Pharma (www.quaypharma.com)

FlavorActiV Limited, Sanderum House, Oakley Road, Chinnor, Oxon, OX39 4TW, UK. Tel No. +44(0) 1844 354 154 email: enquiries@flavoractiv.com; website: www.flavoractiv.com